

FEEDSCAN WEB | FOR CONSULTANTS

GLOSSARY OF TERMS

KEEP IT ALL STRAIGHT WITH THIS QUICK REFERENCE:

Account The name of the account the plot samples should be associated with and billed to. To ensure this is populated correctly, your login should be established and associated with all accounts within your jurisdiction prior to uploading a plot sample spreadsheet.

Additional Emails Include emails of any additional team members who need to receive the results of the plot samples analyzed.

Date Sampled This defaults to the date of sample information generation and is autofilled. This can be adjusted using the calendar function if crafting a large batch of sample labels for a future date of sampling.

'Excluded from Reporting' Notes that are for the laboratory staff only. These will only be seen by the Rock River Laboratory data entry team and are not included in reports.

Feeder The farm, company, or individual associated with the Rock River Laboratory account.

Feeder Farm The farm that is utilizing the feedstuffs to be sampled. This may be the farm associated with the account, or it may be a different farm.

FeedScan Code The random series of numbers and letters assigned to a sample when utilizing the FeedScan app or FeedScan web interface. This code is replaced with a "lab number" when the sample is received by the lab.

Grower Farms The farms growing the feeds to be purchased.

Lab # The unique identifier assigned to a sample upon receipt by the lab. It follows the sample through the analysis process and is included on the final report. The Lab # is different from the FeedScan code.

FEEDSCAN WEB | FOR CONSULTANTS

GLOSSARY OF TERMS

Package The analysis the user would like to run on the sample. This is broken down between the various analysis formats - Near Infrared (NIR), Wet Chemistry, Wet Chemistry - Additional (these must be added to a NIR or Wet Chemistry package), Anti-nutritional Analysis, and Water Analysis.

Rolling Average A feedstuff's analysis results, averaged over time. This is calculated based on the number of samples selected (i.e. "3" will average 3 recent samples from the selected feedstuff).

Sample Descriptions The titles given to each of the feedstuffs that are to be sampled. This description is determined by the user, and may include categorical names, dates, etc.

Submit Multiple Samples Submit all samples for a given day in one fell swoop. Allows users to submit samples in great accuracy - down to the field.

Submit Scheduled Samples If sampling is on a pre-determined schedule, this is outlined within FeedScan Web and users can quickly select from a short list of options of samples that are 'due' for their selected time frame. Sample labels can then be ready for use by the nutritionist or dropped off at the farm for the farm to collect the samples, with all the information already incorporated into FeedScan Web.

Vendors This may include any outside companies or brokers who are utilized in purchasing outside feeds. Users may list up to three vendors in the Sample Submission box.